

# **2003 NMR User Training Course**

National Program for Genomic Medicine High-Field NMR Core Facility,  
The Genomic Research Center, Academia Sinica

09/29-09/30, 2003

**09/30, 2003 Course Handout**

**Software Training: NMRPipe & Sparky**

by

Wen-Jin Wu

## NMRPipe

- **Read the original paper first**

Delaglio, F., Grzesiek, S., Vuister, G.W., Zhu, G., Pfeifer, J. and A. Bax (1995) NMRPipe: A multidimensional spectral processing system based on UNIX pipes. J. Biol. NMR 6, 277-293.

- **Usage**

Processing of 1D-4D NMR data.

- **Availability&Installation**

Free for Academia purpose. NMRPipe website: <http://spin.niddk.nih.gov/bax/software/NMRPipe/>

- **Platform**

Unix, Linux, NT-NMRPipe (with fee, <http://www.resdesigns.com/>)

- **Command Help**

`nmrPipe -fn command -help`,

for example: “`nmrPipe -fn SP -help`” to find out more about the command SP.

Or, simply use “`man SP`”, “`man bruk2pipe`” etc.

## NMRDraw

- **Come with NMRPipe. For viewing NMRPipe processed data. Useful for getting phasing parameters or solvent subtraction scheme for NMRPipe processing script.**

- Some Links for using NMRPipe/NMRDraw:

<http://lysine.pharm.utah.edu/nmrlab/nmrpipe.html>

NMRPipe mailing list:

<http://groups.yahoo.com/group/nmrpipe/>

Linux NMR mailing list

<http://www.nmrfam.wisc.edu/Listarchives/linuxnmr-archive/threads.html>

Note on NMRPipe/Sparky:

[http://www.pharmacy.umaryland.edu/PSC/NMR/proc\\_anal/Pipe2Sparky.html](http://www.pharmacy.umaryland.edu/PSC/NMR/proc_anal/Pipe2Sparky.html)

[http://www.cm.utexas.edu/hoffman/nmrpipe\\_notes.pdf](http://www.cm.utexas.edu/hoffman/nmrpipe_notes.pdf)

Varian data:

[http://www.pharmacy.umaryland.edu/PSC/NMR/proc\\_anal/](http://www.pharmacy.umaryland.edu/PSC/NMR/proc_anal/)

---

- **Environment setup for remote display**

1. Open your x-window program (i.e. Exceed for windows, xhost + for Linux/UNIX, disable firewall)
2. Remote login to the work station
  1. Linux/Unix: `ssh <remote hostname>`
  2. Windows: use Putty or other ssh program
3. Setting DISPLAY Variable to for correct color display

```
> setenv DISPLAY <your hostname> : 0.0 ( csh or tcsh shell )  
( Example: setenv DISPLAY Yarrowca : 0.0  
or  
> export DISPLAY=<your hostname> : 0.0 ( sh or bash shell )
```

---

### Using NMRPipe/NMRDraw

#### 2D Processing: HSQC

1. [winston@GOAT 140]\$ setenv DISPLAY wolf:0.0

snp (source NMRPipe, “snp” is just an alias for “source /usr/nmr/NMRPipe/nmrpipe.cshrc” used in IBMS, check with your system administrator).

2. open a xwindow program such as Exceed, and have the display set up.


3. Change to the directory containing the ser file:

```
cd /goat/data/winston/nmr_data/av600/tepcomp_fhsqc/1
```


(Data acquisition parameters for this file: 64 complex point (R+I=128), DQD in  $^1\text{H}$ , STATES-TPPI acquisition mode in  $^{15}\text{N}$ ).

4. type in “bruker” to bring up the conversion utility window (Figure 1) . (type in “varian” for Varian data).

5. Click on the “Read Parameters” button to update acquisition parameters (figure 1)



**Figure 1**



**Figure 2**

6. Change the parameters that are incorrect. In this case, you should obtain something like Figure 2.
7. "Save Script" then "Execute Script" to do the conversion. "Quite" when is done.

```

goat.nmr - PaTTY
[winston@GOAT 1]$ ls
acqu      audita.txt  fid3.com    hsqc_1.ucsf  pdata      test2.fid
acqu2     cag_par    fid.com     hsqc_pipe    pulseprogram  test.fid
acqu2s    cpdprg3    format.temp hsqc_pipe2   scon       winston.dat
acqu      fid2.com    fq0list     hsqc_pipe3   ser        winston_ph.dat
[winston@GOAT 1]$ more fid3.com
#!/bin/csh


bruk2pipe -in /goat/data/winston/nmr_data/av600/tepcomp_fhsqc/1/ser -bad 0.0
-swap -DMX -decim 48 -dspfvs 12 \
  -xN          2048  -yN          128  \
  -xT          1024  -yT          64   \
  -xMODE       DQD   -yMODE       States-TPPI \
  -xSW         3720.238 -ySW         2311.070 \
  -xOBS        600.135 -yOBS        60.818  \
  -xCAR        4.754  -yCAR        114.264 \
  -xLAB        1H     -yLAB        15N   \
  -ndim        2     -aq2D       States \
  -out /goat/data/winston/nmr_data/av600/tepcomp_fhsqc/1/test.fid -verb -ov

sleep 5
[winston@GOAT 1]$ █

```

(Use “man bruk2pipe” for details on the “bruk2pipe” conversion script). (Use “man var2pipe” for Varian data).

8. Open nmrDraw (type this in), and select the file “test.fid”.
9. Under the *Mouse* button, select *1D Horizontal*, and choose X:1, Y:1 (the 1<sup>st</sup> fid).
10. Under *Proc*, choose Auto-Process 1D.
11. Have the *Phasing* “on”. Phase the spectrum and note down the values of ph0 and ph1.



12. Edit the `hsqc_pipe` macro and also input the above phasing parameters ( $ph0=p0$ ,  $ph1=p1$ ) in the proton dimension.
13. Make the macro executable, and place it under the directory containing the converted file "test.fid": `/goat/data/winston/nmr_data/av600/tepcomp_fhsqc/1`.
14. Execute the `hsqc_pipe` script by typing the name directly.

## Processing Status

[winston@GOAT 1]\$ hsqc\_pipe2

PS 128 of 128

PS 2048 of 2048

---

### The “hsqc\_pipe2” script:

---

```
#!/bin/csh
nmrPipe -in ./test.fid \ #Read in the file
| nmrPipe -fn SP -off 0.5 -end 1 -pow 1 -c 1 \ #SP: Sine bell window function
| nmrPipe -fn ZF -auto \ #ZF: zero filling
| nmrPipe -fn FT -auto \ # FT: Fourier Transformation
| nmrPipe -fn PS -p0 -50 -p1 0.0 -di -verb \ #phase, delete imaginary
#| nmrPipe -fn EXT -left -sw \ #Extract region
| nmrPipe -fn TP \ #Transpose axis (xy to yx)
| nmrPipe -fn SP -off 0.5 -end 1 -pow 1 -c 1 \ #process on the 2nd dimension
| nmrPipe -fn ZF -auto \
| nmrPipe -fn FT -auto \
| nmrPipe -fn PS -p0 0.0 -p1 0.0 -di -verb \
nmrPipe -out ./winston_ph.dat -ov # Write out processed data
```

---

(use “man command Name”, or “nmrPipe -fn Command Name -help” to find out more about each command !)

---



**nmrPipe -fn SP [-off offset] [-end end] [-pow pow]**  
| nmrPipe -fn SP -off 0.5 -end 1 -pow 1 -c 1 \

- Use “nmrPipe -fn SP –help” for help on SP.

**SP: Adjustable Sine Window. [SINE]**

**-off offset [0.0] Sine Start\*PI. (Q1) ;use 0.5 for a Cosine starts at 1 . Sine(0.5 $\pi$ )=1**  
**-end end [1.0] Sine End\*PI. (Q2) ;use 1 for a Cosine ends at 0. Sine( $\pi$ )=0**  
**-pow exp [1.0] Sine Exponent. (Q3) ;1 for Sine, 2 for Sine square.**  
**-c fScale [1] Point 1 Scale.**

- **NMRPipe Processing Functions**

A list of the nmrPipe functions used in the data processing:

**EXT** extracts a region from the current dimension with specified limits  
**FT** applies a real or complex forward or inverse Fourier transform  
**HT** applies a Hilbert transform to reconstruct imaginary data  
**LP** linear prediction  
**POLY** When used in frequency-domain, applies polynomial baseline correction .  
When used in the time-domain, applies solvent correction  
**PS** applies the zero- and first-order phase corrections  
**REV** reverse data order in given dimension  
**SOL** applies solvent correction

**SP** applies sine-bell apodization

**TP/YTP** exchanges data vectors from x and y axis of the data stream

**ZF** applies zero-filling

**ZTP** exchanges data vectors from x and z axis of the data stream

### ***Generic Arguments***

The following arguments are used by more than one function in the shell scripts above.

**-di** deletes imaginary data after the given processing function is performed

**-hdr** extracts parameters recorded during previous processing from spectral header

**-in** specifies the input file or file template

**-inPlace** specifies replacement of the input data by the output result (use with care)

**-inn** activates the inverse mode of a given function


**-out** specifies the output file or file template

**-ov** permits overwriting of any pre-existing files.

**-sw** updates the spectral width and other ppm calibration information


**-verb** permits processing in verbose mode, with status messages

---



**NMRDRAW Display**

*(Since “frequency-jump” was used in this experiment, the 1H-dimension needs to be calibrated.)*



**Sparky Display**

**Use NMRDraw to inspect the spectrum**

- Check phasing of column and row vector by using “1D Horizontal, or 1D vertical” under “Mouse”.
- **Convert to Sparky format (2D):** `pipe2ucsf winston_ph.dat hsqc_1.ucsf`
- Open the spectrum in Sparky:

### 3D Data Processing using NMRPipe/NMRDraw

- change to the directory containing the ser file:  
`/goat/data/winston/nmr_data/nmrfam/hnco_omtky3/1`
- `[winston@GOAT 1]$ brucker`

NMRPipe Conversion Utility Version 97.027.12.56

Spectrometer Input:

Output Template:

Output Script:

Input Protocol:

Output Protocol:

Dimension Count:

Temperature (K):

	x-axis	y-axis	z-axis
Total Points R+I:	<input type="text" value="2048"/>	<input type="text" value="80"/>	<input type="text" value="64"/>
Valid Points:	<input type="text" value="1024"/>	<input type="text" value="40"/>	<input type="text" value="32"/>
Acquisition Mode:	<input type="text" value="DQD"/>	<input type="text" value="States"/>	<input type="text" value="States"/>
Spectral Width Hz:	<input type="text" value="5482.456"/>	<input type="text" value="1729.244"/>	<input type="text" value="1886.437"/>
Observe Freq MHz:	<input type="text" value="500.132"/>	<input type="text" value="125.779"/>	<input type="text" value="50.684"/>
Center Position PPM:	<input type="text" value="4.706"/>	<input type="text" value="175.549"/>	<input type="text" value="118.573"/>
Axis Label:	<input type="text" value="1H"/>	<input type="text" value="13C"/>	<input type="text" value="15N"/>

Read Parameters Save Script Execute Script Hide Script Clear Script Update Script Quit Help

## The ft\_xyz.com nmrpipe macro:

\*\*\*\*\*

```
#!/bin/csh
#
# 3D States-Mode HN-Detected Processing.

xyz2pipe -in fid/test%03d.fid -x -verb \
| nmrPipe -fn SP -off 0.5 -end 0.95 -pow 1 -c 1.0 \
| nmrPipe -fn ZF -auto \
| nmrPipe -fn FT -verb \
| nmrPipe -fn PS -p0 0.0 -p1 0.0 -di \
| nmrPipe -fn EXT -x1 10.0ppm -xn 7.0ppm -sw \
| nmrPipe -fn TP \
| nmrPipe -fn SP -off 0.4 -end 0.95 -pow 1 -c 0.5 \
| nmrPipe -fn ZF -auto \
| nmrPipe -fn FT \
| nmrPipe -fn PS -p0 0.0 -p1 0.0 -di \
#| nmrPipe -fn CS -rs 1.7ppm -sw \
| nmrPipe -fn POLY -auto -ord 0 \
| nmrPipe -fn TP \
| nmrPipe -fn POLY -auto \
| nmrPipe -fn TP \
| pipe2xyz -out ft/B%03d.ft2 -y -ov

xyz2pipe -in ft/B%03d.ft2 -z -verb \
```

```
| nmrPipe -fn SP -off 0.5 -end 0.95 -pow 1 -c 0.5      \  
| nmrPipe -fn ZF -auto                                  \  
| nmrPipe -fn FT -verb                                  \  
| nmrPipe -fn PS -p0 0.0 -p1 0.0 -di                 \  
#| nmrPipe -fn CS -ls 0.2ppm -sw                      \  
| nmrPipe -fn POLY -ord 0 -auto                        \  
| pipe2xyz -out ft/C%03d.ft2 -z -ov
```

\*\*\*\*\*

(The “ft\_xyz.com” should be executable.)

### Processing status message

[winston@GOAT 1]\$ ft\_xyz.com

XYZ2Pipe Partition: Plane 1 to 64 of 64

FT 5104 of 5120

FT 5120 of 5120

XYZ2Pipe Partition: Plane 1 to 128 of 128

FT 71760 of 71808

FT 71808 of 71808

## Convert the Processed Data to Sparky-Readable Format

**3D Data:** For a processed hncO spectrum (3D data) saved as 2D planes with names hncO001.ft, hncO002.ft, ... the command to produce a single 3D NMRPipe file looks like:

```
% xyz2pipe -in hncO%03d.ft -x > hncO_3D.ft
```

you will then use the pipe2ucsf utility to convert the “single file” 3D NMRPipe data to Sparky format file:

```
% pipe2ucsf hncO_3D.ft hncO_3D.ucsf
```

## Sparky

- **Utility:**

View and analyze multidimensional NMR data.

- **Availability & Installation:**

Free @ <http://www.cgl.ucsf.edu/home/sparky/>

- **Platform:**

Unix, Linux, MS Windows on PCs

- **Useful Links:**

(1) Sparky's website

<http://www.cgl.ucsf.edu/home/sparky/>

(2) Sparky Practical:

<http://dolphin.chem.uu.nl/~henry/cursus/sparky.html>

## Using Sparky

### Reading Files:

(Sparky's manual: <http://www.cgl.ucsf.edu/home/sparky/manual/files.html>)

- **Data processed with Bruker XWINNMR:**

To convert Bruker processed data 1/pdata/1/2rr to UCSF format:

```
% bruk2ucsf 1/pdata/1/2rr noe150.ucsf
```



- **Data processed with NMRPipe:**

**(1). 2D Data:** Given NMRPipe data noe150.pipe (2D data) you convert it to UCSF format with:

```
% pipe2ucsf noe150.pipe noe150.ucsf
```

**(2). 3D Data:** For a processed hncO spectrum (3D data) saved as 2D planes (with names hncO001.ft, hncO002.ft, .....hncO0032.ft), first “pack” all files into a single 3D file:

```
% xyz2pipe -in hncO%03d.ft -x > hncO_3D.ft
```

, then use the pipe2ucsf utility to convert the “single” 3D NMRPipe file to Sparky-formatted file:

```
% pipe2ucsf hncO_3D.ft hncO_3D.ucsf
```

### Sparky Software Demonstration

**lbd\_3d\_13\_1 (CBCANH)**

**lbd\_3d\_12\_1 (CBCA(CO)NH)**

Contour level adjustment, changing plan, integrate, vector display, strip plot etc.