

# 2023 NMR Users Training (II)


Basic NMR SOP for  
Small Molecules & Metabolomics Analysis

## Data Collection & Processing using Topspin

Chi-Fon Chang  
Manager, HFNMRC  
Senior Research Specialist/Scientist, GRC  
Academia Sinica  
2023.10.24

# QA1: Do I need high field for better resolution ?

NMR field	<= 500MHz	>= 600 MHz
chemical shift(ppm)	the same	
coupling constant(J, Hz)	the same	
Sensitivity*	lower	higher
Resolution**	lower	better



$$N_{\beta}/N_{\alpha} = \exp(-\Delta E/kT) = \exp[(\gamma\hbar B_0)/(2\pi kT)]$$


$$S/N = \frac{N\gamma_{exc}T_2(\gamma_{det}B_0)^{3/2}\sqrt{ns}}{T}$$

- $S/N$  = signal to noise ratio
- $N$  = number of spins in the system (sample concentration)
- $\gamma_{exc}$  = gyromagnetic ratio of the excited nucleus
- $\gamma_{det}$  = gyromagnetic ratio of the detected nucleus
- $ns$  = number of scans
- $B_0$  = external magnetic field
- $T_2$  = transverse relaxation time (determines the line width)
- $T$  = sample temperature

\* Higher Sensitivity is needed for low concentration sample

\*\* Better resolution is needed for overlap peaks

# 1D <sup>1</sup>H Spectrum in ppm




Higher the field, better the sensitivity and resolution

## QA2: Which Probe could provide information I need?

Probe Type	Regular	Cryoprobe
Sensitivity* (organic solvent)	1	~4
(aqueous solution)	1	~2.5

Probe Coil	Observe (ex: TBO, BBO, Dual)	Inverse (ex: TXI, TCI)
	X-nuclei observed Experiment (1D $^{13}\text{C}$ , 1D $^{31}\text{P}$ )	$^1\text{H}$ -detected Experiments (1D $^1\text{H}$ , 2D COSY/TOCSY 2D HSQC/HMBC )



## NMR & Probes in HFNMRC

NMR & Probe	Topspin	1H (EB)	Others
<b>NEO500_IBMS</b> (Cryo QNP)	<b>TP4.2</b>	2,000	1,000 (13C) 988 (31P) 1,000 (19F)
<b>AVIII600_IBMS</b> (Cryo TCI_005)	TP2.x TP3.x	5,700	710(13C)
<b>NEO600_IBMS</b> (Cryo TCI_121)	<b>TP4.1</b>	6,530	950(13C)
<b>AV600_CHEM</b> (regular TBO)	TP2.x	406	377(13C) 241(31P) 342(19F)
<b>AV800_IBMS</b> (Cryo TXI)	TP2.x	6,200	N/A
<b>AVIII850_IBMS</b> (Cryo TCI)	TP2.x TP3.x	8,500	1,600(13C)

Regular Probe*	1H (EB)	Others
<b>500MHz_TXI</b>	450	
<b>600MHz_TXI</b>	1,218	
<b>600MHz_BBO</b>	465	465 (13C)
<b>600MHz_QXI</b>	1,193	85(31P)
<b>600MHz_TBO</b>		
<b>800MHz_TXI</b>	2,077	

CRYO Probe*	1H (EB)	Others
<b>500MHz_TXI</b>	4,196	
<b>500MHz_QNP</b>	2,000	1,000 (13C) 988 (31P)
<b>600MHz_TCI_005</b>	5,700	710(13C)
<b>600MHz_TCI_121</b>	6,530	950(13C)
<b>800MHz_TXI</b>	6,200	
<b>850MHz_TCI</b>	8,500	1,600(13C)

\*Signal to Noise (S/N) value @ installed date

# Data Collection using Topspin

**Simple Operation Guide for HFNMRC Users**

by Dr. Chi-Fon Chang for small molecules (2022.03.10 updated)


# Topspin Software

	HFNMRC	Version	Upgrade in 2024
1	NEO500	Topspin4.2	Topspin4.3
2	AVIII600	Topspin2.1 Topspin3.2	
3	AV600_CHEM	Topspin2.1	
4	NEO600	Topspin4.1	Topspin4.3
5	AV800	Topspin2.1	Topspin2.1 Topspin3.2
6	AVIII850	Topspin2.1 Topspin3.2	Topspin4.3

	GRC	Version
1	AV600R	Topspin2.1
2	AV600L	Topspin2.1 Topspin3.1

# What's new for Topspin4.x

➤ NEO must use Topspin4.x or higher version




- ◆ Routinely used commends are the same as Topspin2.x or Topspin3.x  
指令與Topspin2.x or Topspin3.x 相通
- ◆ Interface is different but not too difficult to follow  
介面稍有不同, 熟悉即可
- ◆ HFNMR Standard Experiments and SOP are the same  
實驗設定方式與本核心其它NMR相同



# Simple Operation Guide for HFNMRC Users

by Dr. Chi-Fon Chang for small molecules (2022.03.10 updated)

## PART I: Preparation (Software & Hardware)



<Lock> The  $^2\text{H}$  signal is collected by “lock channel” that operates in parallel with the principle channels. “Lock” maintain the center of  $^2\text{H}$  resonance at a constant frequency.


- Deuterated solvents are used to generate the signal to be detected and monitored by the lock system. The frequency and strength of this signal will depend on the solvent used.
- The lock system uses a receiver to monitor this deuterium frequency and makes adjustments to the magnetic field strength accordingly.
- The deuterium frequency is measured several thousand times per second. Hence, as long as the system is locked, the user can be confident that the field is maintained at a constant strength during acquisition.

**\*Important: you must lock the correct solvent, otherwise, the spectrum chemical shift might be incorrect.**

# Simple Operation Guide for HFNMR Users



by Dr. Chi-Fon Chang for small molecules (2022.03.10 updated)

## PART I: Preparation (Software & Hardware)



**<Wobble>** Wobble is to carry out “tuning” and “matching” simultaneously. Type “wobb” for old probe, “atma” or “atmm” for probe with ATM

- Tuning involves adjusting the probe circuitry so that the **frequency** at which it is most sensitive is the relevant transmission frequency (SFO1, SFO2 etc.)
- Matching involves ensuring that the **maximum amount of the power** arriving at the probe base is transmitted up to the coil which lies towards the top of the probe.




**\*Important:** for QNP or BBO probe, you should wobble for the X-nuclei you like to observe

# Simple Operation Guide for HFNMRC Users


by Dr. Chi-Fon Chang for small molecules (2022.03.10 updated)

## PART I: Preparation (Software & Hardware)



**<Shimming>** Shimming is a process in which minor adjustments are made to the magnetic field until the **field homogeneity** (uniformity) is optimized. Improving the homogeneity will result in better spectral resolution.

- Shimming is to adjust the shim coil circuits (ex: Z, Z2, Z3...X,Y....) which are designed to create small magnetic fields that will **cancel out inhomogeneity** in the  $B_0$  magnetic field.




- Type "topshim" for automation, manually shim might be needed

# Simple Operation Guide for HFNMRC Users

by Dr. Chi-Fon Chang for small molecules (2022.03.10 updated)

## PART II: Experiment Set up & Data Collection



<getprosol> getprosol is to read in pulse parameters saved in the “prosol”

- prosol is a file (table) containing Probe/Solvent parameters (calibrated using standard sample, ex: Urea/Methanol ) for the probehead installed on a specific NMR system.

The screenshot shows the 'edprosol' software window. At the top, it says 'Saved Observe and Saved Decouple Prosol Parameter Set for:'. Below this, there are fields for 'Probe:' (Z44896\_0121 CP TCI 600S3 H-C/N-D-05 Z) and 'Solvent:' (generic). There are also dropdown menus for 'Observe' and 'Decouple' channels, both currently set to '1H'. Below these are text boxes for 'Observe Comment:' (Default 1H obs 600) and 'Decouple Comment:' (Default 1H dec 600). At the bottom, there is a table with tabs for '90 deg. Pulses', 'HR Square Pulses', 'HR Shape Pulses', and 'Others'. The table has columns for 'Observe' and 'Decouple', each with sub-columns for 'Nucleus', 'Pulse Width[μs]', 'Att. Lvl.[dB]', and 'Set'. The data in the table is as follows:


Observe				Decouple			
Nucleus	Pulse Width[μs]	Att. Lvl.[dB]	Set	Pulse Width[μs]	Att. Lvl.[dB]	Set	Nucleus
1H	8.00	-7.32		8.00	-7.32		1H
2H	68.00	-14.81		68.00	-14.81		2H
13C	12.00	-19.55		12.00	-19.55		13C
15N	35.00	-19.23		35.00	-19.23		15N

**\*Important: the pulse parameters in “HFNMRC standard parameter set” won’t be the same for current probe. Thus, you must “getprosol” to read in the correct parameters for current probe you are using.**



# Simple Operation Guide for HFNMR Users

by Dr. Chi-Fon Chang for small molecules (2022.03.10 updated)

## PART II: Experiment Set up & Data Collection





`<pulsecal>` is a command to **determine 90 degree pulse** for your sample



- Write down the **power level (dB)** and **corresponding pulse-length(us)**
- Instead of “getprosol” only, you should type  
“getprosol 1H `<90° pulselength>` `< corresponding power level>` “




**\*With optimized 90 degree pulse, the spectrum quality could be better. Especially for multi-dimensional experiments or experiments with specific pulse angles (ex: dept90, dept135, cosy45 )!**

# 90° vs. 30° Pulse



**Signal Intensity=Full** →


**But need to wait longer time for next scan**



**Signal Intensity=1/2** →

**But back to Z faster, can collect more scan**



# Data Processing using Topspin



# NOTE: Topspin file structure 檔案編排方式


Data Folder / Experiment Number / Process Number

Data 檔名 / 實驗編號 (EXPNO) / 圖譜處理編號 (PROCNO)





# Workflow for Data Processing (1)



**SI:  $> 2 * TD$**   
(Size of spectrum)

**SR: (check)**  
(Spectrum Reference)

Use internal or external standard (ex: TPP, DSS) to do calibration

**WDW: EM or QSINE**  
(Window Function)

**1D**

Bruker TopSpin 3.6.0 on NB-000112 as cfchang

Start Process Analyse Publish View Manage

Create Dataset Find Dataset Open Dataset Paste Dataset Read Pars.


ProteinTrain.2019 1 1 G:\data

Spectrum ProcPars AcqPars Title PulseProg Peaks Integrals Sample Structure Plot Fid

S 12 ME

Parameter	Value	Description
SI	32768	Size of real spectrum
SF [MHz]	600.1400000	Spectrometer frequency
OFFSET [ppm]	12.71319	Low field limit of spectrum
SR [Hz]	0	Spectrum reference frequency
HZPPT [Hz]	0.293438	Spectral resolution
SPECTYP	UNDEFINED	Type of spectrum e.g. COSY, HMQC, ...
<b>Window function</b>		
WDW	EM	Window functions for trf, xfb, ...
LB [Hz]	1.00	Line broadening for em
GB	0	Gaussian max. position for gm, 0<GB<1
SSB	2	Sine bell shift SSB (0,1,2,...)
TM1	0	Left limit for tm 0<TM1<1
TM2	0	Right limit for tm 0<TM2<1
<b>Phase correction</b>		
PHC0 [degrees]	98.824	0th order correction for pk

# Workflow for Data Processing (1)



Broker TopSpin 3.6.0 on NB-000112 as cfchang

Start Process Analyse Publish View Manage

Proc. Spectrum Adjust Phase Calib. Axis Pick Peaks Integrate Advanced

Protein Train.2019 2 1 G:\data

2D


Reference	F2	F1	Frequency axis
SI	2048	256	Size of real spec
SF [MHz]	600.1400000	60.8116580	Spectrometer fre
OFFSET [ppm]	10.70628	135.00130	Low field limit of
SR [Hz]	0	0	Spectrum refere
HZpPT [Hz]	3.521259	8.314709	Spectral resoluti
SPECTYP	HSQC		Type of spectrur
Window function			
WDW	QSINE	QSINE	Window function
LB [Hz]	0.30	0.30	Line broadening
GB	0	0.1	Gaussian max. p
SSB	2	2	Sine bell shift SS
TM1	0	0.1	Left limit for tm 0
TM2	0	0.9	Right limit for tm

Protein Train.2019 11 2 G:\data [3]


3D

Reference	F3	F2	F1
SI	2048	32	128
SF [MHz]	600.1400000	60.8116580	150.9053230
OFFSET [ppm]	10.70628	135.00130	70.00701
SR [Hz]	0	0	0
HZpPT [Hz]	1.595570	66.517670	37.741547
AQORDER	3-2-1		
SPECTYP	UNDEFINED		
Window function			
WDW	QSINE	QSINE	QSINE
LB [Hz]	1.00	0.30	0
GB	0.1	0.1	0
SSB	2	2	2
TM1	0.1	0.1	0.1
TM2	0.0	0.0	0.0


# Workflow for Data Processing (2)




**1D : efp**  
(Fourier Transfer on 1 dimension)




**2D : xfb**  
**3D: xfb & ft3d**  
(Fourier Transfer on "both", ie, 2 dimension)



# Workflow for Data Processing (3)



1D → .ph  
(PHC)/PHC1




Bad Phase !!



Good Phase !!




# Workflow for Data Processing (3)




2D → select 2 or more “rows” → adjust PHC0 & PHC1 similar to 1D

**Step 1:** pick peaks by click right mouse and add, then click “row” to phase the row



**Step 2:** use mouse to adjust phase 0 (zero order) and phase 1 (first order)




Save the good phase | + - ||| ≡ ⚙

**Step 3 :**Exit the “phase mode”




# Workflow for Data Processing (3)



2D → select 2 or more “rows” → adjust PHC0 & PHC1 similar to 1D



Bad Phase !!



Good Phase !!



# Workflow for Data Processing (4)




## Other useful Tips

**STSR / STSI : adjust window**  
(adjust spectrum window)

**ME\_mod: (ex: LPfr)**  
(Linear prediction)  
**NCOEF: 8\*n (ex: 8, 16...)**

**REVERSE: FALSE / TURE**  
(check from spectrum)



# Hands-on

## Let's try Some Experiments

- 1D 1H one pulse (zg)
- 1D 1H solvent suppression (noesypr1d)
- 1D 13C DEPT (deptq)
- 2D 1H-1H COSY
- 2D 1H-13C HSQC