


# 2004 NMR Users Training Course

## *02/27 Lab Session Work Sheet*

### Simple Outline

- Hardware (CRYOPROBE)  
Please make sure everything is normal for the Cryoprobe system !!
- Loading Sample
  - (i) check temperature ( “edte” {command line}: edit temperature )



( Temp limitation for CryoProbe is 283K-333K, check Temp Cal. Table for setting desired temp )

- (ii) check spinner (sample holder with yellow line) position
  - (iii) lift gas (make sure there's lift gas) {click button}
  - (iv) load sample
  - (v) lift off (now the sample should go down slowly) {click button}
  - (vi) lock the solvent ( “lock” {command line}, and click on your D-solvent)
- Wobble ( Don't need to change cable anymore!!! )
    - (i) check “edasp” setting {command line}
    - (ii) type “wobb” {command line} ( this might take 20-30 sec, please be patient)
    - (iii) type “acqu” {command line} to observe wobble curve
    - (iv) adjust “Tune” and “Match” of probe ( **Carefully!! Never over tune !!**)
    - (v) click on “ wobb-sw” to switch to another nuclei if necessary ( F3→F2→F1)
  - Shimming (gradient shim)

## Exercises

- Data Collection Easy 1-2-3
  1. start a new experiment : [new](#) or [edc](#)
  2. read in parameters for standard experiment : `rpar "CRP_*`
  3. Adjust parameters for your own sample(p.s.) Users just need to optimize experiments by adjusting parameters that list under "AQ parameter to check"
- Pulse & Parameter Calibration 1-2-3
  1. optimize O1 , especially for H2O suppression purpose
  2. optimize useful 1H pulse length  
( [Please read the Pulse and Power Recommendation from Bruker](#) )
  3. `d1+aq` has to be [longer than 1 sec](#)

## Experiment (1)

Experiment Name: [2D 15N-1H HSQC](#)

Experiment Type: [Using echo-antiecho, f1: H, f3:N](#)

Standard Parameter Set: [CRP\\_2D\\_15N\\_HSQC\\_ETSI](#)

Pulse Program: [hsqcetf3gpsi2](#)

AQ parameters to check

1H pulses

[p11](#)(high power, ex: -5dB), [p1](#)(90° pulse at p11), [p28](#) (trim pulse, ex:1m)

Others

[cnst4](#) ( $J_{H-N}$ , ex: 90Hz), [d24](#) ( $1/4J_{H-N}$ ), [o1](#) (for 1H), [o3](#) (for 15N),

[1 sw](#), [1 td](#) (for F1 dimension, N), , [2 sw](#), [2 td](#) (for F2 dimension, H)

[d1](#), [ns](#)(=1\*n), [ds](#)(>=16), [rg](#)

---

## **\*\* Let's see what happen without 15N decouple??**

Experiment Name: [2D 15N-1H HSQC](#)

Experiment Type: [Using echo-antiecho, f1: H, f3:N](#)

Standard Parameter Set: [CRP\\_2D\\_15N\\_HSQC\\_ETSI](#)

Pulse Program: [hsqcetf3gpsi2](#)

AQ parameters to check

Keep everything the same as above, but [pL16=120dB](#) !!!!!

## Experiment (2)

Experiment Name: 2D 15N-1H TROSY

Experiment Type: Using echo-antiecho, f1: H, f3:N

Standard Parameter Set: CRP\_2D\_15N\_TROSY\_ETSI

Pulse Program: trosyef3gpsi

AQ parameters to check

1H pulses

p11(high power, ex: -5dB), p1(90° pulse at p11)

Others

cnst4 ( $J_{H-N}$ , ex: 90Hz), o1 (for 1H), o3 (for 15N),

1 sw, 1td (for F1 dimension, N), , 2 sw, 2 td (for F2 dimension, H)

d1, ns(=4\*n), ds(>=16), rg

---

## Experiment (3)

Experiment Name: 2D 15N-1H CRINEPT

Experiment Type: Using echo-antiecho, f1: H, f3:N

Standard Parameter Set: CRP\_2D\_15N\_CRINEPT

Pulse Program: crineptgpph

AQ parameters to check

1H pulses

p11(high power, ex: -5dB), p1(90° pulse at p11)

p11 ( shape pulse , ex: 1000us), sp1 ( shape pulse level for Sinc1.1000)

Others

cnst4 ( $J_{H-N}$ , ex: 90Hz), o1 (for 1H), o3 (for 15N),

1 sw, 1td (for F1 dimension, N), , 2 sw, 2 td (for F2 dimension, H)

d1, ns(=2\*n), ds(>=16), rg

(p.s.) Using AU-program “ splitterinept” to process data

---